

Targeting Vector Preparation for Electroporation

- 1) Plasmid DNA should be prepared using EndoFree Plasmid Maxi protocol.
- 2) Digest 100 ug of targeting vector with a restriction digestion enzyme that cuts the vector outside arms of homology overnight (take ~1/50 of sample for gel analysis later).
- 3) Do ethanol precipitation for a couple of hours with 2.5 volume of 100% ethanol on ice (or -20 degree).
- 4) Centrifuge at maximum speed for 10 minutes to pellet DNA.
- 5) Wash DNA with 70% ethanol. Centrifuge at maximum speed for 5 minutes. Repeat wash step and allow pellet to air-dry.
- 6) Resuspend pellet in sterile PBS to a concentration of 1ug/ul.
- 7) Run a 0.8% gel with 500 ng samples of uncut, digested & unpurified, and digested & purified along with 1ug **λ -HindIII Digest Marker** (NEB cat. #N3012) to check for linearization and verify DNA concentration.

IMPORTANT! Project will not be accepted without λ -HindIII Digest Marker gel. Nanodrop readings can be misleading. Comparison to the marker allows the most accurate measurement of vector concentration. This will maximize the targeting efficiency.

- 8) Provide gel picture and 50ug of linearized DNA for electroporation.